Applying laser technology for cancer care Health Tuesday 2021: Cancer innovations



Lasers and optics for personalized medicine and better life.


Seppo Orsila 2nd Feb 2021

Modulight, Inc.

- The only vertically integrated medical laser manufacturer in the world
- Profitable privately held company with > 20 years track record
- Exclusive supplier to 10+ pharmas & world leading cancer centers in USA & Europe
- Long track record of making all kinds of semiconductor lasers & pumps (400–2000 nm)
- Certified ISO13485 QMS in 2011, ISO9001 & 14001 in 2002
- Co-branded & OEM products currently shipping to all continents except Antarctica
- 600 000 pcs of 9xx-nm pulsed lasers in 2015 for LIDAR



Photoimmunotherapy (PIT) concept



PDT: combination treatment in oncology



Photoimmunotherapy

- Less invasive targeted treatment of tumours
- Reducing and inhibiting side effects
- Treatment of tumour margins in surgery
- Initiating or boosting anticancer immunoresponse







Examples of photoimmunotherapy application areas

- □ Head & neck cancers
- □ Esophageal cancers
- □ Bile duct cancer
- Eye melanoma
- Non-small lung cell carcinoma
- Mesothelioma
- Glioblastoma
- Prostate
- □ Many others...

Topics today

Photoimmunotherapy in oncology Laser therapy for glioblastoma treatment (German, Japan) Laser therapy for uveal melanoma (US) Laser therapy for lung cancer (Tampere)



Modulight: Cure for GBM


- Modulight has been working with the best experts in the world since 2015 to create a cure for glioblastoma and solve traditional challenges associated with laserbased treatments
- Initial results of glioblastoma iPDT are very promising and show
 5-year survival increase to over
 50% with iPDT (<10% with standard of care treatment methods)
- Modulight has developed a state-ofthe art laser platform with cloud connection and real-time treatment monitoring capabilities for GBM phototherapy

Photodynamic therapy for glioblastoma

- 5-Aminolevulinic Acid (5-ALA) accumulates in tumor and undergoes a convertion to protoporphyrin IX (PpIX).
- PpIX has an absorption peak at 635 nm, and in the presence of appropriate light source, its fluorescence acts with bimodal function:
 - **5**-ALA-PDT

 \rightarrow series of irreversible photochemical and photobiological events that cause directly damage and killing tumor cells

- □ Fluorescence diagnostic marker
- → assists neurosurgeons to visualize the extent and margins of tumors
- → Photosensitizer monitoring before, during and after treatment



PDT treatment scenarios for GBM





Photodynamic therapy


treatment

Potential recurrence block via immunoresponse



modulight

Photoimmunotherapy



Benefits of novel Modulight laser

- Same laser system and same fibers used for illumination and detection
- Addition of treatment monitoring is made effortless by integrating them to the treatment flow
- Eight channels providing dosimetry data on light distribution over tumor area
- Photosensitizer presence in tissue is confirmed prior to starting the treatment
- Integrated sensors & cloud-based analytics that help in real-time to guide the surgery
- Collected diagnostic data can be viewed from <u>cloud.modulight.com</u> during the treatment and downloaded later



Topics today

0

Photoimmunotherapy in oncology Laser therapy for glioblastoma treatment (German, Japan) Laser therapy for uveal melanoma (US) Laser therapy for lung cancer (Tampere)



Uveal melanoma is cancer of the eye

- Kills aggressively: 5-year survival rate only 15%, half develop metastases, 5 % of all melanomas are uveal
- Traditional treatment is heavy: even the removal of the eye

Modulight ML7710 laser with slit lamp adapter

- Mountable to Zeiss SL120/SL130/SL115C and Haag Streit type slit lamps
- Set up, calibration and training by Modulight field engineers
- Direct field and service support from the factory
- Multispot treatment supported with adjustable spot sizes
- Easy to use touch screen graphical user interface

The spot diameter can be adjusted from 0.5 mm to up to 20 mm with a contact lens. It is a range beyond compare.

Dr. Carol Shields Wills Eye Hospital

Investigational Drug Trial

Topics today

0

Photoimmunotherapy in oncology Laser therapy for glioblastoma treatment (German, Japan) Laser therapy for uveal melanoma (US) Laser therapy for lung cancer (Tampere)



Photodynamic therapy for intrabronchial tumors

- ML7710 operated at 665 nm and sodium chlorin E6 photosensitizer
- Illumination dose for the target area is 50–75 J/cm²
- The aim in most patients has been a palliative recanalization of the bronchial lumen to alleviate symptoms like breathlessness and hemoptysis

PDT of intrabronchial obstructive tumors

- Before the first PDT treatment the patient's main bronchus was almost completely blocked by the tumor
- Approximately one week after the treatment bronchus is partially open and the collapse of the lung is prevented





- Second PDT treatment was completed a few months after the first one.
- Pictures on the left show PDT treatment area before and a few weeks after the second treatment



Summary

- PDT/photoimmunotherapy is a proven modality for cancer treatment; latest advances include treatment monitoring before, during and after treatment, targeted drugs, combination treatments, and immunotherapeutic response
- Clinical trials for glioblastoma are expected to give results in a few years
- Photodynamic therapy my be game changer to uveal melanoma as a first line of treatment
- Intrabronchial tumor PDT treatment expected to save more lives in the future and improve the level of palliative treatment

Thank you! modulight

MORE GREAT STORIES BY FOLLOWING US IN TWITTER, FACEBOOK & YOUTUBE:

MORE STORIES TWITTER: @ORSILA @MODULIGHT









noduligh

